



STEINFURTH®
Electromechanical Measuring Systems



Product Range for Quality Control

New beverage types and ingredients, packaging materials and designs, closure solutions, price and cost sensibility are the sources of the actually dynamic changing keys to the consumer satisfaction.

The beverage and packaging quality can't be seen or monitored these days anymore separately – both areas of those are the inseparable parts of the QA influencing each other.

Steinfurth Core Competencies:

Strong focus on solutions for combined quality control on beverages and beverage packages.

Most optimal product performances by cooperation between our customers and our engineering team.

Cost and time efficiency as result of innovation and leadership in engineering services for beverages and packaging industry.



Steinfurth CPA

Compact Package Analyzer as customized extendable, modular Mini Lab Concept for combined quality control on beverages and beverage packages (CO₂ content, torque, fill height, brix, pH, conductivity, ...)



Steinfurth CO₂MS/CDA

Carbon Dioxide Analyzer for automatic CO₂ measurement and CO₂ shelf life studies on packaged beverages



Product Range for Quality Control

Product Range for Quality Control



Steinfurth CDA OST

Carbon Dioxide Analyzer for non-destructive CO₂ measurement in the area of CO₂ Shelf Life Studies on packaged beverages (PET & Glass bottles)



Steinfurth TMS 2000/TorqTraQ

Torque Measurement System for manual and mobile torque measurement on bottle closures and capping heads



Steinfurth CCT/TCT

Cup/Tank Carbonation Tester for standardized, handheld CO₂ measurement on beverage dispensing systems, tanks and KEGs



Steinfurth DV

Digi Vacu for handheld vacuum and pressure measurement on closures, bottles & cans



Steinfurth BPT/CPT

Bottle/Can Pressure Tester for handheld pressure measurement on PET-bottles, glass bottles with crown cork closures and cans



Steinfurth TS92.5

Pressure Calibrator for mobile and flexible calibration of manometers and various pressure measuring instruments




Steinfurth TMS 5010/4010/4000

Torque Measurement System for user independent, automatic torque measurement on bottle closures



Steinfurth Silver Bottle TS CDA-Silverbottle

Silver Bottle for standardized checking of pressure and temperature calibration on CO₂ testers based on Henry-Dalton physical law



Product Range for Quality Control

Product Range for Quality Control



Steinfurth Orange Bottle

Orange Bottle for standardized checking of torque calibration on automatic and manual torque testers and operator training on manual torque testers



Steinfurth Rotech KEG

Logger system for monitoring of the KEG washing/sterilization and filling process in KEG filling lines




Steinfurth SF-PastControl

Logger system for monitoring of the pasteurization process in tunnel and chamber pasteurizers



Steinfurth LCS 710P

Laboratory Carbonation System for single sample carbonation studies and preparation of tasting samples in bottles and cans



Steinfurth SF-Control B System

Pressure and Temperature Logger system for monitoring of pressure and temperature in containers



Steinfurth SFT

Foam Tester for automatic foam stability measurement on beer foam



Steinfurth SF-CleanControl

Logger system for monitoring of the cleaning process on returnable bottles in bottle washer machinery



Steinfurth Sampler

Sampling device for CO2 loss & safety protected beverage sample extraction out of glass bottles, PET-bottles and cans



Product Range for Quality Control

Product Range for Quality Control



Steinfurth AS380C

Automatic Sampler for multiple CO₂ loss & safety beverage sample extraction out of glass bottles and cans



Steinfurth SMART

Side Wall Abuse Resistance Tester for automated & standardized side wall impact test on beverage cans and metal bottles



Steinfurth AS380CP

Automatic Sampler for multiple CO₂ loss & safety beverage sample extraction out of PET-bottles



Steinfurth 2 CDA-MK6

Carbon Dioxide Analyzer for automatic CO₂ measurement and CO₂ shelf life studies on packaged beverages - twin container version



Steinfurth LabShaker

Laboratory Shaker/Tumbler for beverage sample preparation and homogenization



Steinfurth OPT 180

Opening Performance Tester for standardized, automatic opening performance test on bottle closures and finishes

Visit Us Online:
www.steinfurth.com

Your Steinfurth Team



Steinfurth - Headquarters

Steinfurth Mess-Systeme GmbH

Elektromechanische Mess-Systeme
Bonifaciusring 15
45309 Essen
Germany

Phone +49-201-85517-0

Telefax +49-201-85517-20

Internet: <http://www.steinfurth.de>

E-Mail: info@steinfurth.de

Steinfurth, Inc.

Electromechanical Measuring Systems 11390
Old Roswell Rd.
Ste. 126
Alpharetta, GA 30009
USA

Phone +1-678-500-9014 & +1-404-918-5061

Telefax +1-678-840-7744

Internet: <http://www.steinfurth.com>

E-Mail: info@steinfurth.com

